

INTERSPACE® KNEE SPACER

TEMPORARY KNEE SPACER WITH GENTAMICIN

INDICATIONS FOR USE:

InterSpace is indicated for temporary use (maximum 180 days) as a total knee replacement (TKR) in skeletally mature patients undergoing a two-stage procedure due to a septic process.

InterSpace is applied on the femoral condyles and on the tibial plate following removal of the existing implant and radical debridement. The device is intended for use in conjunction with systemic antimicrobial antibiotic therapy (standard treatment approach to an infection).

InterSpace is not intended for use for more than 180 days, at which time it must be explanted and a permanent device implanted or another appropriate treatment performed (e.g. resection arthroplasty, fusion etc.). Because of the inherent mechanical limitations of the device material (gentamicin/polymethylmethacrylate), InterSpace Knee is only indicated for patients who will consistently use traditional mobility assist devices (e.g. crutches, walkers, canes) throughout the implantation period.

PRODUCT DESCRIPTION:

InterSpace is composed of two articulating independent elements: a tibial component composed of a flat base upon which a femoral component articulates. Both components must be fixed to the bone with bone cement, such as Cemex Genta. InterSpace is composed of fully formed gentamicin/polymethylmethacrylate (PMMA) radiopaque bone cement. The two piece design mimics an ultra-congruent condylar knee prosthesis. The design allows for motion of the leg for basic mobility (sitting, standing, walking) under limited weight bearing conditions (e.g. crutches, walkers, canes). InterSpace is made with bone cement containing 1.87% gentamicin base.


CONTRAINDICATIONS:

Use of InterSpace is contraindicated in the following situations:

- The patient's condition is such that a two-stage arthroplasty procedure is contraindicated due to decreased immune response or other relevant systemic clinical conditions.
- Bone loss precluding adequate support of the prosthesis.
- Lack of adequate competence (anatomical and functional) of peripheral ligamentous apparatus and extensor mechanism.
- The procedure is unjustified due to deficiencies in the patient's muscular, nervous or vascular systems.
- Poor bone quality (as in osteoporosis) could cause the prosthesis to migrate or to fracture host bone.
- Infection of the TKR cannot be confirmed.
- The infected TKR devices cannot be removed.
- The infecting pathogens are resistant to gentamicin.
- The patient is sensitive (allergic) to gentamicin, aminoglycosides or PMMA bone cement.
- A systemic or secondary remote infection is suspected or confirmed.
- The patient does not have a TKR and the infection is secondary to trauma, septic arthritis or other surgical procedures.
- The patient does not have sufficient bone stock to allow insertion and fixation of the prosthesis
- The patient has neuromuscular disorders that do not allow control of the knee joint.
- The patient's weight, age or activity level would cause the surgeon to expect early failure of the system.
- The patient is unwilling or unable to use protected weight bearing mobility throughout the implantation period (e.g. crutches, canes, walkers etc).
- Myasthenia gravis

MICROBIOLOGY:

Gentamicin sulphate is an aminoglycoside antibiotic derived from the actinomycete *Micromonospora purpurea*. The molecular weight is 449.55. The product contains no preservative or sodium bisulfite. Gentamicin sulphate is a complex of the gentamicins C1, C1a and C2, which illustrated below:



Gentamicin	Mol. Formula	R1	R2	R3
C1	C ₂₁ H ₄₃ N ₅ O ₇	CH ₃	CH ₃	H
C1a	C ₁₉ H ₃₉ N ₅ O ₇	H	H	H
C2	C ₂₀ H ₄₁ N ₅ O ₇	H	CH ₃	H
C2a	C ₂₀ H ₄₁ N ₅ O ₇	H	H	CH ₃
C2b	C ₂₀ H ₄₁ N ₅ O ₇	CH ₃	H	H

INFORMATION REGARDING GENTAMICIN SULPHATE RELEASED FROM PMMA:

Mechanism of action

Bacteria tend to adhere to surfaces where they can multiply and create a defensive barrier called a biofilm (complex structure mainly made by extracellular polysaccharides and proteins). Bacteria embedded in a biofilm are more resistant to most antibiotic therapy because the glycoprotein structure is difficult for antimicrobial agents to penetrate. PMMA, due to its surface characteristics, is one of the materials with the highest risk of bacterial colonization. It has been demonstrated *in vitro* that the presence of antibiotics in PMMA reduces bacterial adhesion.

Gentamicin activity is primarily directed against aerobic, gram negative bacilli. The action against most gram positive bacteria is limited. Gentamicin is active against sensitive strains of enterococci and streptococci at concentrations which can be achieved clinically only when combined with a penicillin. Gentamicin is active *in vitro* against more than 90% of strains of *S. aureus* and 75% of *S. epidermidis*. Gentamicin has been shown to be active against most strains of the following organisms both *in vitro* and in clinical infections.

Most common susceptible pathogens

Gram positive bacteria

Staphylococcus aureus; *Streptococcus pyogenes*; *Streptococcus pneumoniae*;
Streptococcus faecalis; *Listeria monocytogenes*

Gram negative bacteria

Citrobacter; *Enterobacter*; *Escherichia coli* *Klebsiella* spp.; *Proteus mirabilis*; *Proteus vulgaris*; *Morganella morganii*; *Providencia* spp.; *Salmonella* spp.; *Serratia*; *Shigella* spp.; *Pseudomonas aeruginosa*

Bibliography

Goodman & Gilman's The Pharmacological Basis of Therapeutics 12th Edition, March 2011, Chapter 54 (McDougall C, Chambers) - McGraw Hill, New York

PHARMACOLOGICAL WARNINGS:

In vitro elution studies (microbiological method) has shown that the amount of gentamicin released never exceeds the recommended systemic adult dose for gentamicin (5-7 mg/kg/day) according to the Goodman and Gilman's recommendations (adults with normal renal function). It is therefore unlikely that the amount of gentamicin released from the device and absorbed locally will result in serum levels in the toxic range.

Nonetheless trough concentrations which exceed 2µg/ml for longer than 10 days have been associated with toxicity (systemic administration). InterSpace should be used with caution (mainly in the first day of implantation of the spacer) in conjunction with other nephrotoxic or ototoxic drugs. The use of antibiotic-loaded bone cement for fixation of the InterSpace device may increase the potential of toxic drug reactions. The InterSpace device should be used with caution in patients who are predisposed to or who have preexisting clinical conditions that would put them at risk for gentamicin

toxicity (e.g. renal dysfunction, dehydration, advanced age etc.). All patients should be monitored for toxic blood levels of gentamicin, nephrotoxicity and ototoxicity while the device is in situ. This is especially important in elderly subjects and in those receiving other nephrotoxic and/or ototoxic drugs.

POTENTIAL ADVERSE EVENTS:

The following serious and frequent adverse effects may be associated with use of the InterSpace device. Although some effects are not directly attributable to the device itself, the surgeon should be aware of these potential complications and be ready to treat the patient accordingly.

General Surgical Risks

- venous thrombosis
- transitory hypotension
- myocardial infarction
- pulmonary embolism
- arrhythmias
- sudden death

TKR Surgery Risks

- damage to femur or tibia
- damage to blood vessels
- nerve damage, bone bed damage
- arthrofibrosis
- limb length discrepancy
- phlebitis, thrombophlebitis
- hematoma
- wound healing problems
- extensive blood loss

InterSpace Device Risks

- gentamicin toxicity:
ototoxicity; nephrotoxicity
- PMMA sensitivity
- recurrent infection
- device breakage
- device dislocation
- device loosening
- debris release
- foreign body reaction
- difficulty in removing the device

ATTENTION: Since the device must be fixed with antibiotic bone cement, the surgeon must be aware of its negative effects. Recurrences of infections, although rare, have been known to recur even with IV antibiotic use. All aminoglycosides have the potential to produce reversible and irreversible vestibular, cochlear and renal toxicity. Adverse reactions to Gentamicin Sulphate are not expected at the low levels used within InterSpace. However the following adverse reactions have been associated with larger doses, typical of prescribed dosages of Gentamicin Sulphate for systemic parenteral administration.

High serum peaks of aminoglycoside caused by once-daily drug administration are well tolerated; the once-daily regimens are just as safe as or safer than multiple-dose regimens.

Neurotoxicity

- Manifested as both auditory and vestibular ototoxicity, including irreversible hearing loss
- Numbness
- Skin tingling
- Muscle twitching
- Convulsions

Neurotoxicity - Adverse effects on both the vestibular and auditory branches of the eighth nerve have been noted, especially in patients receiving high doses or prolonged therapy, in those given previous courses of therapy with an ototoxic drug, and those suffering from dehydration.

Symptoms include dizziness, vertigo, tinnitus, roaring the ears and hearing loss. Hearing loss is usually irreversible and is manifested initially by diminution of high-tone acuity. Gentamicin and tobramycin closely parallel each other in regard to ototoxic potential.

Nephrotoxicity

- Usually in patients with pre-existing renal damage
- Also in patients with normal renal function to whom aminoglycosides and administered for longer periods or in higher doses than recommended
- The symptoms of which may manifest after cessation of therapy

Nephrotoxicity - Renal function changes, as shown by rising BUN, NPN, and serum creatinine and by oliguria, cylindruria, and increased proteinuria, have been reported, especially in patients with a history of renal impairment who are treated for longer periods or with higher doses than those recommended. Adverse renal effects can occur in patients with initially normal renal function.

Clinical studies and studies in experimental animals have been conducted to compare the nephrotoxic potential of gentamicin and tobramycin. In some of the clinical studies and in the animal studies, tobramycin caused nephrotoxicity significantly less frequently than gentamicin.

In some other clinical studies, no significant difference in the incidence of nephrotoxicity between tobramycin and gentamicin was found. Neuromuscular blockage or respiratory paralysis, more commonly in patients with myasthenia gravis or Parkinson's Disease. In general aminoglycosides have little allergenic potential; both anaphylaxis and rash are unusual. Rare hypersensitivity reaction - including skin-rashes, eosinophilia, fever, etc. – have been reported. Other reported adverse events possibly related to gentamicin include: anemia/granulocytopenia, thrombocytopenia, fever, rash, exfoliative dermatitis, itching, urticaria, nausea, vomiting, diarrhea, headache, lethargy, mental confusion and disorientation. Laboratory abnormalities possibly related to gentamicin include increased serum transaminases including AST and ALT, increased serum LDH and bilirubin, decreased serum calcium, magnesium, sodium and potassium; and leukopenia, leukocytosis, and eosinophilia.

PATIENT PRECAUTIONS:

The physician must instruct the patient as follows:

- Protected weight bearing mobility must be used throughout the implantation period (e.g. crutches, canes, walkers etc).
- Any condition that tends to impose severe loading on the InterSpace device should be avoided (e.g. obesity, participation in active sports, unprotected weight bearing, likelihood of falls etc.).
- Report any pain, discomfort or trauma with the affected limb.
- The InterSpace product must be explanted after the temporary use.

Because the InterSpace device was designed for temporary implantation under protected load bearing conditions, the patient should be periodically evaluated with X-ray with respect to the condition of the interface fixation between bone and device, bone trophism and other relevant clinical conditions during the rehabilitation phase.

USE DURING PREGNANCY AND BREAST-FEEDING:

There are no tests that demonstrate the utilization safety of InterSpace during pregnancy, breast-feeding. Knee revision surgery should be avoided during the first three months of pregnancy. This product is indicated for applications in the remaining gestation period only when it is believed impossible to save the joint or preserve the patient's life through other forms of intervention.

USE IN CHILDREN:

There are no tests that demonstrate the InterSpace is safe to use in children. The device should only be used in skeletally mature individuals.

PRECAUTIONS FOR USE:

Familiarity with and attention to appropriate surgical techniques for knee arthroplasty revision surgery and familiarity with proper use of the InterSpace device is essential for successful use of the device. Only surgeons who have reviewed the surgical technique regarding InterSpace implantation and are aware of the limitations of its application should utilize the device.

The user must not modify the device in any way, including not adding other antibiotics as the effects pharmacologically and structurally cannot be predicted. The user must not allow damage to the device. Any alteration or damage to the component may reduce fatigue strength and may result in failure under load. The wear rate of prosthesis component contact surfaces is greatly accelerated if loose fragments of bone, bone cement, or other particulate debris become detached and act as an abrasive in the articular and modular interfaces.

The expected useful life of the InterSpace component may be compromised in a very large or overweight individual and/or one who does not adequately protect the amount of activity and weight placed on the knee. It is recommended to always implant the largest component size possible. It is essential that the patient use mobility assist devices (e.g. crutches, walker) during the implantation period.

Aggressive assembly methods are not required for proper implantation of the device. During the application do not subject the device to excessive forces (e.g. hammer strikes) that could cause damage. Any damage to the device may affect the fatigue strength and lead to failure under load. A fine needle aspiration from the joint site and antibiotic susceptibility testing should be performed prior to implantation of InterSpace. All patients should be instructed on the limitations of the prosthesis and the need for a subsequent surgery to implant a definitive prosthesis. Patients should be taught to govern their activities accordingly. Post-operative care is important. Implants must not be reused. Any implant, once used, should be discarded even though it may appear undamaged. Failure to adhere to these recommenda-

tions will result in increased probability of poor function, loosening, wear, fracture or premature failure. Do not use the InterSpace device in cases where the existing implant components cannot be completely explanted. Do not use the InterSpace implant in joints that contain osteosynthesis devices that could mechanically interfere with its function. Do not allow the component to remain implanted for more than 180 days. The device was tested to be safely used for not more than 6 months. If this period is extended for too long this can lead to wearing, development of debris and eventually to breakage that can cause pain, inflammation and bone reabsorption. After removal of the InterSpace device, the wound site should be thoroughly irrigated to remove all bone cement debris prior to implantation of a definitive prosthesis or alternative surgical procedures (e.g. resection arthroplasty, fusion etc.). Failure to remove cement and/or bone debris may shorten the survival of the revision implant. Using InterSpace under conditions other than the indicated use is unlikely to provide benefit to the patient and increases the risk of the development of drug-resistant bacteria.

Warning: Do not re-sterilize and/or re-use the device. It is designed for single-use on a single patient. Re-sterilization should not be carried out since it can cause infection risks for the patient. Re-sterilization can also alter the morphology, the efficiency of the antibiotics and the mechanical features of the device, causing a malfunction of the same with serious risks for the patient's health. Re-using the device after its extraction must by all means be avoided since it can cause contamination and worsening of the patient's infection. During the extraction, the spacer can also be damaged or residues of cement can be left on the device.

UTILIZATION, IMPLANTATION AND EXPLANTATION:

It is important to maintain strictly aseptic surgical techniques. Avoid washing with aqueous solutions the spacer before or after implantation to maintain optimal levels of antibiotic release.

Selection of the proper size InterSpace depends on the judgement of the surgeon with relationship to the requirements of the patient. The surgeon shall become thoroughly familiar with the technique of implantation of the prostheses by: (1) reading appropriate literature, (2) training in the operative skills and techniques required for InterSpace Knee arthroplasty revision surgery, and (3) reviewing information regarding use of instrumentation for sizing, implantation and explantation of the component.

For the selection of the size, transparent radiograph overlays and InterSpace Trial devices are available.

Caution: the thickness generated by the cement used for the fixation is not reflected in these templates.

In order to prevent dislocation, the same measures taken for a permanent total knee replacement (TKR) are advised, plus other specifics such as:

- 1) Choice of the correct size (see below);
- 2) Proper cement fixation (Cemex Genta) of the components;
- 3) Insertion with appropriate joint tension of the soft tissues around the knee joint;
- 4) Application of an orthopaedic abduction brace, possibly articulated, to assist with flexion for cases with greater risk of dislocation.

INTERSPACE SIZE CHOICE:

The most suitable size is selected on the basis of:

- a) Dimensions of the removed implant;
- b) The type of bone defect.
- c) Ligamentous apparatus state;
- d) Flexion and extension gaps.

Among the sizes, the size to be implanted is that which:

- is nearest to the size of the removed implant; and
- achieves the best compromise between stability and joint mobility during the operation. It is important to achieve stable full extension and 90° flexion.

Once the size has been chosen, it is advisable to apply both components without cement and reduce the joint to evaluate stability and joint function. The two components must be fixed with Cemex Genta bone cement filling the gaps of the bone-device interface. The entire contact surface of the components must be cemented to create continuity between spacer and bone. Missing or insufficient cement - especially in the posterior portion of the condyles - may weaken the structure of the device.

Apply the femoral part first, and wait for the polymerization of the cement, then proceed with the application of the tibial part. Ensure that excess cement does not adhere to the articulating surfaces.

Reduce the joint before the cement of the tibial part has cured to allow for the self-centering of the tibial component in relation to the femoral one. Keep the joint in full extension until complete polymerization has occurred to allow for the tibial component fixation. Clean the area of any debris.

When the suture and the reconstruction of the extensor apparatus has been completed the knee must be stable, but not too tight, and have a joint extension ranging from 0° to 90°.

Post-operative treatment

As a general rule, post-operative treatment is superimposable with a primary knee prosthesis, with the difference that the weight-bearing can be only partial (use of crutches).

Partial weight-bearing must be assessed on an individual basis in relation to the anatomic conditions of the femur and tibia, bone trophism and the clinical conditions of the patient during rehabilitation stages.

In particular, excessive weight-bearing and forced mobilization should be avoided to minimize the risk of InterSpace damaging the biological structure.

If the surgeon deems it necessary, on the basis of the stability achieved and, especially, the condition of the extensor apparatus, temporary use of an articulated post-operative orthosis may be prescribed.

Explantation

The InterSpace device is not intended for use as a permanent prosthesis and must be removed within 180 days of implantation. Osteotomes, mallets and other revision instruments may be used to aide in the explantation procedure. Care should be taken to assure that the wound site is thoroughly cleaned of all bone cement debris prior to implantation of a definitive prosthesis or performing an

alternative surgical procedure (e.g. resection arthroplasty, fusion etc.). Failure to remove cement and/or bone debris may shorten the survival of the revision implant.

DISPOSAL


Disposal of the device should be in accordance with local waste regulations.

HOW SUPPLIED:

InterSpace implants are supplied sterile. Do not re-sterilize. Prior to use, all packages should be inspected for integrity. If a package is damaged, opened or contaminated in any way, it must not be used.

CAUTION:

Federal law restricts this device to sale by or on the order of a physician.



REF	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)	Gentamicin Base
SPK0022	54	40	60	36	16	0.9 g
SPK0122	64	47	70	42	17	1.3 g
SPK0222	74	54	80	48	18	1.8 g
SPK0322	84	61	90	54	19	2.7 g

Symbols:



USE BY



STERILIZED USING ETHYLENE OXIDE



DO NOT RESTERILISE



BATCH
NUMBER



CONTENTS STERILE UNLESS PACKAGE
IS DAMAGED OR OPENED



CATALOGUE
NUMBER



DO NOT RE-USE



CONSULT INSTRUCTIONS FOR USE



CAUTION



6800 Poplar Avenue | Suite 120 | Memphis, TN 38138

901.453.3141 | info@OsteoRemedies.com | OsteoRemedies.com

OsteoRemedies and the corporate mark are registered trademarks of OsteoRemedies, LLC

InterSpace® is manufactured by Tecres® S.p.A.

Via A. Doria, 6 - 37066 Sommacampagna - Verona, ITALY
and is distributed in North America by OsteoRemedies®

©2023

OST006 Rev 00 10-2023