



# KNEE SPECIFICATIONS



REMEDY®  
FEMORAL COMPONENT



REMEDY®  
TIBIAL  
COMPONENT

REMEDY®  
TIBIAL  
WEDGE/INSERT

## REMEDY® KNEE SPACER

Description	Catalog #	(MM)									Gentamicin Base (g)
		1	2	3	4	5	6	7	8	9	
REMEDY® Tibial Component 60mm	RKTBSM				60	36	25	14	7.8		0.4
REMEDY® Tibial Component 70mm	RKTBMD				70	42	25	14	8.2		0.6
REMEDY® Tibial Component 80mm	RKTBLG				80	48	25	14	8.8		0.9
REMEDY® Femoral Component 54mm	RKFMSM	54	41.6	9.5							0.5
REMEDY® Femoral Component 64mm	RKFMMD	64	49.3	10.5							0.8
REMEDY® Femoral Component 74mm	RKFMLG	74	56.3	11.5							1.2
REMEDY® Femoral Component 84mm	RKFMXL	84	63.4	12.5							1.6
REMEDY® Tibial Wedge/Insert 60mm	RKINSM				60	36		14.5		10	0.3
REMEDY® Tibial Wedge/Insert 70mm	RKINMD				70	42		14.5		10	0.5
REMEDY® Tibial Wedge/Insert 80mm	RKINLG				80	48		14.5		10	0.7