



HIP SPECIFICATIONS


REMEDY® AND REMEDY SPECTRUM™ GV MODULAR HEAD



REMEDY® ACETABULAR CUP




REMEDY®
Acetabular Cup
is indicated for
use with the
REMEDY® Hip
Spacer System.



REMEDY® AND REMEDY SPECTRUM™ GV MODULAR STEM



REMEDY® AND REMEDY SPECTRUM™ GV MODULAR LONG STEM



REMEDY® AND REMEDY SPECTRUM™ GV HIP SPACER

Hip Component Description	REMEDY® Catalog #	REMEDY SPECTRUM™ GV Catalog #	(MM)								REMEDY® Gentamicin Base (g)	REMEDY SPECTRUM™ GV	
			1	2	3	4	5	6	7	8		Gentamicin Base (g)	Vancomycin Base (g)
REMEDY® ACETABULAR CUP 46mm ID/54mm	RHACSM	—	58.5	31							0.3	—	—
REMEDY® & REMEDY SPECTRUM™ GV Femoral Head - 46 mm	RHHDSM	GVHDSM	46	42.3							0.9	0.9	0.9
REMEDY® & REMEDY SPECTRUM™ GV Femoral Head - 54 mm	RHHDMD	GVHDMD	54	50.9							1.6	1.6	1.6
REMEDY® & REMEDY SPECTRUM™ GV Femoral Head - 60 mm	RHHDLG	GVHDLG	60	57.3							2.3	2.3	2.3
REMEDY® & REMEDY SPECTRUM™ GV Femoral Stem - Small	RHSTSM	GVSTSM			111	10	16.5	11.3	35.6	17	0.5	0.5	0.5
REMEDY® & REMEDY SPECTRUM™ GV Femoral Stem - Medium	RHSTMD	GVSTMD			112	11	21.7	15.5	35.6	17	0.6	0.6	0.6
REMEDY® & REMEDY SPECTRUM™ GV Femoral Stem - Large	RHSTLG	GVSTLG			117	11.5	24	16.5	35.6	17	0.7	0.7	0.7
REMEDY® & REMEDY SPECTRUM™ GV Femoral Long Stem - Small	RHLSSM	GVLSSM			227	10	16.5	11.3	35.6	17	0.6	0.6	0.6
REMEDY® & REMEDY SPECTRUM™ GV Femoral Long Stem - Medium	RHLSMD	GVLSMD			227	11	21.7	15.5	35.6	17	0.8	0.8	0.8
REMEDY® & REMEDY SPECTRUM™ GV Femoral Long Stem - Large	RHLSLG	GVLSLG			231	11.5	24	16.5	35.6	17	0.9	0.9	0.9