

REMEDY® & REMEDY SPECTRUM® GV HIP SPACERS**1**

Overview

The REMEDY® Hip Spacer and the REMEDY SPECTRUM® GV Hip Spacer are part of the treatment foreseen in a two-stage procedure performed in the event of permanent prosthesis infection. The REMEDY® Hip Spacer implant and the REMEDY SPECTRUM® GV Hip Spacer are intended for temporary use only (180 days or less). The spacer allows basic joint mobility and releases antibiotic into the joint area to protect the implant from bacterial colonization. A second surgery will be required at a later date to remove the REMEDY® Hip Spacer or the REMEDY SPECTRUM® GV Hip Spacer and replace it with a permanent hip joint implant.



REMEDY® and
REMEDY SPECTRUM®
Modular Heads

REMEDY® and
REMEDY SPECTRUM®
Modular Long Stem

REMEDY® and
REMEDY SPECTRUM®
Modular Stem



Questions

For further product information, please contact Customer Service.

Caution

Federal law restricts this device to sale by or on the order of a physician.

For internal use



Manufactured By: _____



6800 Poplar Avenue | Suite 120 | Memphis, TN 38138

901.453.3141 | info@OsteoRemedies.com | OsteoRemedies.com

OsteoRemedies, REMEDY, SPECTRUM and the corporate mark are trademarks of OsteoRemedies, LLC

©2019 • 0419 V5