




HIP SPECIFICATIONS


REMEDY® AND
REMEDY SPECTRUM® GV
MODULAR HEAD




REMEDY®
ACETABULAR
CUP



REMEDY® AND
REMEDY SPECTRUM® GV
MODULAR STEM



REMEDY® AND
REMEDY SPECTRUM® GV
MODULAR LONG STEM



REMEDY® AND REMEDY SPECTRUM® GV HIP SPACER

Hip Component Description	REMEDY® Catalog #	REMEDY SPECTRUM® GV Catalog #	(MM)								REMEDY®	REMEDY SPECTRUM® GV		
			1	2	3	4	5	6	7	8	Gentamicin Base (g)	Gentamicin Base (g)	Vancomycin Base (g)	
REMEDY® ACETABULAR CUP 46mm ID/54mm OD	RHACSM	—	58.5	31	54							0.3	—	—
REMEDY® & REMEDY SPECTRUM® GV Femoral Head - 46 mm	RHHDSM	GVHDSM	46	42.3								0.9	0.9	0.9
REMEDY® & REMEDY SPECTRUM® GV Femoral Head - 54 mm	RHHDMD	GVHDMD	54	50.9								1.6	1.6	1.6
REMEDY® & REMEDY SPECTRUM® GV Femoral Head - 60 mm	RHHDLG	GVHDLG	60	57.3								2.3	2.3	2.3
REMEDY® & REMEDY SPECTRUM® GV Femoral Stem - Small	RHSTSM	GVSTSM			111	10	16.5	11.3	35.6	17		0.5	0.5	0.5
REMEDY® & REMEDY SPECTRUM® GV Femoral Stem - Medium	RHSTMD	GVSTMD			112	11	21.7	15.5	35.6	17		0.6	0.6	0.6
REMEDY® & REMEDY SPECTRUM® GV Femoral Stem - Large	RHSTLG	GVSTLG			117	11.5	24	16.5	35.6	17		0.7	0.7	0.7
REMEDY® & REMEDY SPECTRUM® GV Femoral Long Stem - Small	RHLSSM	GVLSSM			227	10	16.5	11.3	35.6	17		0.6	0.6	0.6
REMEDY® & REMEDY SPECTRUM® GV Femoral Long Stem - Medium	RHLSDM	GVLSMD			227	11	21.7	15.5	35.6	17		0.8	0.8	0.8
REMEDY® & REMEDY SPECTRUM® GV Femoral Long Stem - Large	RHLSLG	GVLSLG			231	11.5	24	16.5	35.6	17		0.9	0.9	0.9