

# REMEDY<sup>®</sup>

## ACETABULAR CUP SPACER



### COMPLEX PROBLEM SIMPLE SOLUTION

- The first available pre-formed PMMA, antibiotic eluting acetabular cup spacer
- 54mm OD cup articulates with 46mm head to improve device stability and mitigate dislocations
- Helps protect acetabular integrity during staged procedures
- Reduce OR time vs. intra-op molded spacers
- REMEDY<sup>®</sup> Hip System offers complete modularity
  - 54mm OD / 46mm ID cup
  - 3 head sizes
  - 6 stem sizes
  - 17mm of head translation

OsteoRemedies, LLC | 6800 Poplar Avenue | Suite 120 | Memphis, TN 38138

1-800-OSTEO-XL | 901-453-3141 | info@OsteoRemedies.com | OsteoRemedies.com



REMEDY, REMEDY SPECTRUM, OsteoRemedies and the corporate mark are trademarks of OsteoRemedies, LLC • REMEDY and SPECTRUM are manufactured by OsteoRemedies, LLC  
©2019 • MKT018 • 11-19-V1