

REMEDY® ACETABULAR CUP

1

Overview

The REMEDY® & REMEDY SPECTRUM® GV Hip Spacer are temporary hip implants that are used as part of a two-stage procedure. If you are receiving the REMEDY® or REMEDY SPECTRUM® GV Hip Spacer, your surgeon has recommended a two-stage procedure for you because of an infected hip implant. In the first stage, the REMEDY® or REMEDY SPECTRUM® GV Hip Spacer will be surgically implanted in your hip. The REMEDY® and REMEDY SPECTRUM® GV Hip Spacer are intended for temporary use only (180 days or less). It allows basic joint mobility and releases antibiotic into the joint area to protect the implant from bacterial colonization.

