

R E S O R B A B L E B E A D K I T (R B K)

OSTEOBOOST[®] SELECT[™]

100% PURE, SYNTHETIC,
FULLY RESORBABLE BONE VOID FILLER
WITH NANO-CRYSTALLINE TECHNOLOGY

OVER 20,000 CLINICAL CASES



Highly Moldable
Putty Packs Into Voids
and can be Used to Coat
Implant Surfaces

**Bioengineered Triphasic Graft,
Enhanced with nHA**



- Calcium Sulfate
- Beta-tricalcium pPhosphate
- Nano-crystalline Hydroxyapatite (nHA)

Biocompatible

Neutral pH

Optimized Resorption Rate

Indicated for
Gaps and Voids Located
in Extremities, Spine,
Pelvis, or Cranium



Complete resorption in 42 days¹

Beads can be
Packed Into Gaps and
Dead Space

OSTEOBOOST® SELECT™ RESORBABLE BEAD KIT

Complex Problem - Simple Solution

OSTEOBOOST® SELECT™ RBK utilizes a time-tested bone void filler with over 10 years clinical use to fill voids and gaps not intrinsic to the stability of the bone structure.



Sterile Resorbable Bead Convenience Kit

- Bone Void Filler Powder
- Distilled Water
- Mixing Spatula
- Spreading Spatula
- Bead Mat

Elongated Bead Shape Maximizes Surface Area



X-Small
(3 x 4.5mm)



Small
(4.5 x 7mm)




Medium
(5.5 x 9.5mm)

OSTEOBOOST® SELECT™ RBK's unique elongated bead shape is available in three bead sizes and one double-sided bead mat for both 5cc and 10cc kits.

Front of Mat with Small and Medium Beads

Back of Mat with X-Small Beads



OsteoRemedies, LLC | 6800 Poplar Avenue | Suite 120 | Memphis, TN 38138

1-800-OSTEO-XL | 901-453-3141 | info@OsteoRemedies.com | OsteoRemedies.com

OsteoRemedies, OsteoBoost and the corporate mark are registered trademarks of OsteoRemedies, LLC
1. Data on file. ©2024 • MKT032 Rev 04 02-2024



OSTEOREMEDIES®
ADVANCED MEDICAL TECHNOLOGIES