

R E S O R B A B L E B E A D K I T (R B K)

OSTEOBOOST[®] SELECT[™]

100% PURE, SYNTHETIC,
RESORBABLE BONE
VOID FILLER

WITH NANO-CRYSTALLINE TECHNOLOGY



Highly Moldable
Putty Packs Into Voids
and can be Used to Coat
Implant Surfaces


**Bioengineered Triphasic Graft,
Enhanced with nHA**

- Calcium Sulfate
- Beta-tricalcium Phosphate
- Nano-crystalline Hydroxyapatite (nHA)

Biocompatible

Neutral pH


Optimized Resorption Rate



OVER 26,000
CLINICAL CASES


Indicated for
Gaps and Voids Located
in Extremities, Spine,
Pelvis, or Cranium

Months to Resorption



Optimized resorption rate with
dissolution in 42 days¹

Beads can be
Packed Into Gaps and
Dead Space



OSTEOBOOST® SELECT™ RESORBABLE BEAD KIT

Complex Problem - Simple Solution

OSTEOBOOST® SELECT™ RBK utilizes a time-tested bone void filler with over 10 years clinical use to fill voids and gaps not intrinsic to the stability of the bone structure.

Sterile Resorbable Bead Convenience Kit

- Bone Void Filler Powder
- Distilled Water
- Mixing Spatula
- Spreading Spatula
- Bead Mat



Elongated Bead Shape Maximizes Surface Area



X-Small
(3 x 4.5mm)



Small
(4.5 x 7mm)



Medium
(5.5 x 9.5mm)

OSTEOBOOST® SELECT™ RBK's unique elongated bead shape is available in three bead sizes and one double-sided bead mat for both 5cc and 10cc kits.

Front of Mat with Small and Medium Beads

Back of Mat with X-Small Beads

Medium
(5.5 x 9.5mm)

Small
(4.5 x 7mm)

X-Small
(3 x 4.5mm)

OsteoRemedies, LLC | 6800 Poplar Avenue | Suite 120 | Memphis, TN 38138

1-800-OSTEO-XL | 901-453-3141 | info@OsteoRemedies.com | OsteoRemedies.com

OsteoRemedies, OsteoBoost and the corporate mark are registered trademarks of OsteoRemedies, LLC
1. Data on file. ©2024 • MKT032 Rev 07 06-2024



OSTEOREMEDIES®
ADVANCED MEDICAL TECHNOLOGIES