

BIOWASH™

S U R G I C A L I R R I G A T I O N +


- ▶ NO RINSE
- ▶ NO EXTRA TIME
- ▶ HIGHEST EFFICACY OF XBIO® TECHNOLOGY

FOR COMPLEX,
HIGH-RISK AND
REVISION SURGERY



BIOWASH™ SURGICAL IRRIGATION+

1. Remove foil seal completely by peeling tab
2. Insert the spike firmly into either round port



POWERED BY: XBIO™



x Only



VIEW IFU

OSTEOREMEDIES
ADVANCED MEDICAL TECHNOLOGIES

Distributed By:
OsteoRemedies, LLC
6800 Poplar Ave, Suite 120 • Germantown, TN
901-453-3141

REF QTY
BW0500 1 (500 mL)

CONTAINS
Sodium Citrate, Citric Acid,
Sodium Lauryl Sulfate, Water



NMS-80146 REV A

INFO@OSTEOREMEDIES.COM

OSTEOREMEDIES.COM



OSTEOREMEDIES
ADVANCED MEDICAL TECHNOLOGIES


BIOWASH™ SURGICAL IRRIGATION⁺

BIOWASH™ Product Features¹

The Only No-Rinse Solution Designed Specifically For High-Risk and Revision Procedures

- Safe for extended exposure
- No discomfort, no irritation
- Suitable for a wide range of patients
- Simple application, no immediate rinsing needed, one step process

Testing of BIOWASH™ was conducted against pathogens which commonly cause orthopedic infections. Compared to irrigants requiring immediate rinsing, BIOWASH™ can safely remain in the surgical wound. BIOWASH™ delivers a cumulative reduction in bioburden contamination, surpassing the effectiveness of the antiseptic irrigants shown below.



In Vitro Testing: Time-to-Kill testing per the duration of intended use.

Catalog #	Description	Quantity	Size
BW0500-04	BIOWASH™	4 Each	500ml Each

INGREDIENTS

Citric Acid, Sodium Lauryl Sulfate, Sodium Citrate

REFERENCES

1. Documents: NMS-40027 Rev A
2. Document: 1167.4-US-en issue date - 2023 - 11
3. Bashyal RK, Mathew M, Bowen E, James GA, Stulberg SD. A Novel Irrigant to Eliminate Planktonic Bacteria and Eradicate Biofilm Superstructure With Persistent

Effect During Total Hip Arthroplasty. J Arthroplasty. 2022 Jul; 37(7S): S647-S652. doi: 10.1016/j.arth.2022.01.045. Epub 2022 Feb 23. PMID: 35210150

NOTE: *In Vitro* results have not been shown to definitively correlate with clinical outcomes.

OsteoRemedies, LLC | 6800 Poplar Avenue | Suite 120 | Memphis, TN 38138
1-800-OSTEO-XL | 901-453-3141 | info@OsteoRemedies.com | OsteoRemedies.com

