

# XPERIENCE<sup>®</sup>

## A HIGHER STANDARD IN IRRIGATION

- » Clinically evaluated
- » No rinse
- » Improved physician workflow
- » Biocompatible
- » Manual or powered irrigation



OSTEOREMEDI<sup>®</sup>  
ADVANCED MEDICAL TECHNOLOGIES



# POWERFUL IRRIGATION BACKED BY SCIENCE



**Cleansing the wound reduces the microbial load and limits the risks.**<sup>1</sup> XPERIENCE solution is designed to rinse away debris and microorganisms from the wound bed and has been featured in multiple published retrospective analyses.



**NO RINSE REQUIRED**



**BIOCOMPATIBLE**



**REDUCES MICROBIAL LOAD**

## Work Flow

Unlike traditional irrigants, XPERIENCE is biocompatible and does not require rinsing. This improves physician workflow by reducing the number of steps and allowing for use throughout the procedure.

## Efficacy

XPERIENCE irrigation solubilizes and removes wound debris including microorganisms.

### Citric Acid

#### *Buffer*

Combines with sodium citrate to maintain an optimal pH.

### Sodium Citrate

#### *Buffer*

Buffers solution to maintain the appropriate pH.

### Sodium Lauryl Sulfate

#### *Surfactant*

Reduces the surface tension within the solution.

1. Natsuhara KM et. al. Mortality During Total Hip Periprosthetic Joint Infection. J Arthroplasty. 2019 Jul;34(7S):S337-S342. doi: 10.1016/j.arth.2018.12.024. Epub 2018 Dec 24.