


REMEDY[®] STEMMED KNEE SPACER**1**

Overview


The REMEDY[®] Stemmed Knee Spacer is part of the treatment in a two-stage procedure performed if you have a permanent prosthesis infection. The REMEDY[®] Stemmed Knee Spacer implant is intended for temporary use only (180 days or less). It allows basic joint mobility and releases antibiotic into the joint area and adjacent to the device to protect the implant from bacterial colonization. A second surgery will be required at a later date to remove the REMEDY[®] Stemmed Knee Spacer and replace it with a permanent knee joint implant.



REMEDY[®]
Femoral Component



REMEDY[®]
Tibial Component



REMEDY[®]
Tibial Insert Wedge



REMEDY[®]
Stem Component