

XPERIENCE[®]

A HIGHER STANDARD IN IRRIGATION

- » Clinically evaluated
- » No rinse
- » Improved physician workflow
- » Biocompatible
- » Manual or powered irrigation



OSTEOREMEDI[®]
ADVANCED MEDICAL TECHNOLOGIES



POWERFUL IRRIGATION BACKED BY SCIENCE



Cleansing the wound reduces the microbial load and limits the risks.¹ XPERIENCE solution is designed to rinse away debris and microorganisms from the wound bed and has been featured in multiple published retrospective analyses.



NO RINSE REQUIRED



BIOCOMPATIBLE



REDUCES MICROBIAL LOAD

Work Flow

Unlike traditional irrigants, XPERIENCE is biocompatible and does not require rinsing. This improves physician workflow by reducing the number of steps and allowing for use throughout the procedure.

Efficacy

XPERIENCE irrigation solubilizes and removes wound debris including microorganisms.

Citric Acid

Buffer

Combines with sodium citrate to maintain an optimal pH.

Sodium Citrate

Buffer

Buffers solution to maintain the appropriate pH.

Sodium Lauryl Sulfate

Surfactant

Reduces the surface tension within the solution.

1. Natsuhara KM et. al. Mortality During Total Hip Periprosthetic Joint Infection. J Arthroplasty. 2019 Jul;34(7S):S337-S342. doi: 10.1016/j.arth.2018.12.024. Epub 2018 Dec 24.